

# BREAKAWAY CHIROPRACTIC TABLE

**medsurge**   
*technics for life*

**CONFY**CARE



## **BREAKAWAY CHIROPRACTIC TABLE**

### **Description**

This is Breakaway Chiropractic table a versatile table providing flexibility for your practice.

Chiropractic drop table use the high-quality durable and easy to clean medical professional PVC and high-density sponge, provide more comfortable environmental health.



### Standard Specification

- One electric motor, standard supports "Timotion" brand for Hi-Lo function.
- Height control from 48.5cm to 66cm in approx 30 seconds by foot switch.
- Three drop sections: Cervical, Thoracic and Lumbar section.
- Abdominal breakaway design: chest section can be raised or lowered, making the treatment approach more diversified.
- Adjustable Headrest featuring single gas struts, angle adjustable from  $-15^{\circ}$  to  $+20^{\circ}$ .
- Adjustable Thoracic section featuring single gas struts, angle adjustable from  $-18^{\circ}$  to  $+20^{\circ}$
- Adjustable Pelvic section featuring manual strut, angle adjustable from  $0^{\circ}$  to  $+20^{\circ}$
- Adjustable foot section featuring manual strut, angle adjustable from  $0^{\circ}$  to  $+50^{\circ}$
- Chiropractic drop table able to lift up to 250kgs working weight from its lowest height.
- Retractable castors 50mm for excellent mobility and safety.
- Chiropractic drop table sleek design with two rounded frame legs and two fixed wheels.
- Adjustable feet on two frame legs work for uneven surfaces.
- Fire/stain/mildew/oil/water/wear resistant Med PVC available in 6 kinds of colors.
- Certificates: CE, FDA, ISO13485-2016, EN60601



- Chiropractic drop table use professional 5cm diameter medical pressure wheel, so that the mobile device more convenient.
- Can be for the waist, neck, shoulders, foot and other parts of the role of adjuvant therapy to do the syndrome.
- Legs unique 2-adjustment screw design, so that the bed for the ability to adapt to improve the ground, more stable.
- Unique tilt, inverted position and lower waist rotation, making the treatment approach more diversified.
- Professional Chiropractic drop table, the use of Dun pressure principle, so that the safety of treatment is guaranteed.

### **Dimensions**

- Upholstery Length: 183cm
- Upholstery Width: 51cm
- Open: 183 x 51 x (48.5-66)cm
- Foam Thickness: 6cm high density foam + 1cm memory soft foam
- Foam Density: 100 g/m<sup>2</sup>
- Packing: 192 x 62 x 51cm (Heavy carton)
- N.W./G.W.: 92/117kgs
- 20GP/40HQ: 32/72pcs





**MEDSURGE HEALTHCARE LIMITED**  
MUSEUM HILL CENTRE, MUSEUM ROAD  
P.O. BOX 75534 00200  
NAIROBI, KENYA  
Tel: 0720 714 337  
Email: [info@medsurgehealth.co.ke](mailto:info@medsurgehealth.co.ke)  
Web: [www.medsurgehealth.co.ke](http://www.medsurgehealth.co.ke)