

ELECTRICAL WHEELCHAIR





-An SGS certified company.







Electrical Wheelchair

Description

- High load bearing aluminum chair frame
- Two type of rear wheel is optional: Steel wire rear wheel / Magnesium alloy five fork wheel
- Drive mode:drive rear wheel
- With electromagnetic brake, more safety for user
- Battery:Lead-acid Battery (2pcs)
- Detatchable armrest & footrest
- Speed: 1-6km/h (Adjustable)













MUSEUM HILL CENTRE, MUSEUM ROAD P.O. BOX 75534 00200 NAIROBI, KENYA Tel: 0720 714 337

Tel: 0720 714 337 Email: info@medsurgehealth.co.ke Web: www.medsurgehealth.co.ke

