


Endotracheal Tubes

medsurge 
technics for life



Endotracheal tubes different sizes

Indication:

Endotracheal tube inserted into the trachea to prevent airway obstruction and respiratory failure.

Features:

- Made of thermoplastic PVC
- X-ray marker over the entire tube length
- Equipped with side murphy eye
- Graduated every 1 cm
- Double depth marker to control the position of the tube
- Low-pressure, high-volume cuff
- Tube size marked in 4 places (on the connector, on the pilot balloon and in two places on the tube body)
- Detachable connector with an outer diameter of 15 mm
- Pilot balloon with lot and size indication
- Latex-free
- Single-use
- EO sterilized
- Individual packaging: paper-foil



MUSEUM HILL CENTRE, MUSEUM ROAD
P.O. BOX 75534 00200
NAIROBI, KENYA
Tel: 0720 714 337
Email: info@medsurgehealth.co.ke
Web: www.medsurgehealth.co.ke