

EYE WASH STATION

medsurge 
technics for life



EYE WASH STATION

Description

The pH Neutral solution will neutralize the chemical immediately to minimize serious injury to the eye - essential for all workplaces with acids, alkali or chemicals. Unit will achieve a 15+ minute rinse with the patented Plum duo eyepiece - 40% of incidents involve both eyes, be protected with the Plum duo eyepiece

Comes with:

- (3) 1000ml sterile Saline bottles with the duo eyepiece
- (1) 200ml pH Neutral eyecup sterile bottle
- (5) 20ml saline ampoules
- (2) sterile eyepads
- (6) wipes

Features

- All contained wall mount cases with mirrors.
- Solution has a 3 year expiration.

Total rinse time:

- c 17 1/2 min. 15 min. (both eyes)
- c. 2 1/2 min. (1 eye) with the pH eyecup and 20ml ampoules

Benefits of Plum First Aid Eye Wash

- Effective, optimal eye rinsing solution: for general eye rinsing or neutralisation in case of contact of acids or alkalis with eyes or body
- Flexibility: a wide range of products to provide the optimal solution for your needs, including stations, wall-mounted solutions in various sizes, and smaller First Aid Eye Washes that can be mounted in a car, put in a tool box or even carried in a belt pouch.
- Sterile and hygienic: when sealed, they have a 3.5-year warranty. The bottles require no maintenance or cleaning. Expiry date is on the label. Easy and quick to use: in case of accident, rinsing of eyes or body part is possible within seconds. Clear and precise instructions for use are on the label.
- High level of functionality: the ergonomic design of the cap fits snugly into the eye socket, ensuring efficient, gentle eye rinsing. Plum First Aid Eye Washes are sterile and CE certified according to the Medical Device Regulations.

Plum eye wash solutions

Your eyesight is one of the most important senses. So take care of your eyes. Regardless how careful you are accidents can happen, but everybody can make sure that they are prepared. With our eye wash solutions, you have done your best.

If an accident happens, seconds can decide whether your eyesight is damaged or not. Therefore, quick access to eye wash is crucial for all workplaces where there is a risk of eye injuries. Our wide product range ensures the best possible solution specifically for your needs. Choose small and large wall-mounted stations directly at the workplace. Or small, mobile solutions for the van, the tool box and the first aid box, and on the cleaning trolley in the special belt bag.

The green Plum Eye Wash deals with dirt. The blue pH Neutral is crucial in connection with accidents involving acids and alkali.

For chemical accidents

If acid or lye comes into contact with the eye, damage to the cornea begins immediately. For such accidents, we recommend our pH Neutral product, which quickly and gently neutralises the liquid in the eye until the normal pH of 7.4% is reached. Further information and documentation on this process is available.

When dirt enters the eye

Time is also of the essence in accidents where metal, wood, dust or any kind of dirt gets into the eye. In such situations, the immediate use of first aid eye wash can prevent the foreign material from sticking to the eyeball, thus avoiding further serious injury. Plum First Aid Eye Wash contains a sterile sodium chloride solution with a pH (0.9%) equal to the pH of the eye fluid. Using Plum First Aid Eye Wash therefore avoids the unpleasant stinging sensation when flushing the eyes and prevents further contamination.



Plum pH Neutral DUO

The bottle contains 4.9% sterile phosphate solution, which quickly neutralizes even concentrated acidic or alkaline substances in the event of contact with the eyes. It can be up to 50 times more effective than pure water. The liquid contains phosphates, which are already present in the human body. The hydrogen phosphate component of the buffer neutralizes the acids, while the other part of the buffer, dihydrogen-phosphate, neutralizes the alkalis. It comes with an easy-to-use bottle, a shower head, a dust cover and clear instructions for use. The DUO design allows two eyes to be washed simultaneously. It is also suitable for stand-alone use or for filling wall-mounted stations. This small, easy-to-handle bottle also fits in a jacket pocket, first aid box, tool bag or a special belt bag. Shelf life: 3.5 years.

The product is available in the following sizes:

- 500 ml bottle (art. no. 4801). Rinse time: approx. 2 minutes. Carton: 6 x 500 ml
- 500 ml bottle (art. no. 4604). Rinse time: approx. 5 minutes. Carton: 12 x 500 ml



Plum pH Neutral

The bottle contains 4.9% sterile phosphate solution, which quickly neutralizes even concentrated acidic or alkaline substances in the event of contact with the eyes. It can be up to 50 times more effective than pure water. The liquid contains phosphates, which are already present in the human body. The hydrogen phosphate component of the buffer neutralizes the acids, while the other part of the buffer, dihydrogen phosphate, neutralizes the alkalis. It comes with an easy-to-use bottle, a shower head, a dust cover and clear instructions for use. It is also suitable for stand-alone use or for filling wall-mounted stations. This small, easy-to-handle bottle also fits in a jacket pocket, first aid box, tool bag or a special belt bag. Shelf life: 3.5 years.

The product is available in the following sizes:

- 200 ml bottle (art. no. 4752). Rinse time: approx. 2 minutes. Carton: 10 x 200 ml
- 1000 ml bottle (art. no. 4746). Rinse time: approx. 5 minutes. Carton: 12 x 500 ml



Plum DUO Eyewash

The bottle contains 0.9% sterile sodium chloride solution. The DUO design for simultaneous washing of both eyes. The Plum first aid eye wash is suitable for rinsing dirt (such as dust, dirt, metal and wood splinters, etc.) and even certain chemicals such as oils and solvents. It comes in an easy-to-use bottle, with an ergonomic cap, a dust cover and clear instructions for use. It is also suitable for stand-alone use or for filling wall-mounted stations. This small, easy-to-handle bottle also fits in a jacket pocket, first aid box, tool bag or a special belt bag. Shelf life: 3.5 years.

The product is available in the following sizes:

- 500 ml bottle (art. no. 4861). Rinse time: approx. 2 minutes. Carton: 6 x 500 ml
- 1000 ml bottle (art. no. 4800). Rinse time: approx. 5 minutes. Carton: 6 x 1000 ml



Plum pH Neutral DUO

The bottle contains 4.9% sterile phosphate solution, which quickly neutralizes even concentrated acidic or alkaline substances in the event of contact with the eyes. It can be up to 50 times more effective than pure water. The liquid contains phosphates, which are already present in the human body. The hydrogen phosphate component of the buffer neutralizes the acids, while the other part of the buffer, dihydrogen phosphate, neutralizes the alkalis. It comes with an easy-to-use bottle, a shower head, a dust cover and clear instructions for use. The DUO design allows two eyes to be washed simultaneously. It is also suitable for stand-alone use or for filling wall-mounted stations. This small, easy-to-handle bottle also fits in a jacket pocket, first aid box, tool bag or a special belt bag. Shelf life: 3.5 years.

The product is available in the following sizes:

- 500 ml bottle (art. no. 4801). Rinse time: approx. 2 minutes. Carton: 6 x 500 ml



Wound and Eyewash spray

With Plum Wound and Eyewash Spray you get immediate care to minor accidents such as cuts and scrapes on the skin or dust and dirt in the eyes. It helps you to minimize the risk of further infection or damage. It is a portable, easy to use, saline-based and pain-free spray. It works at any angle, providing a balanced flow rate of sterile 0.9 % saline solution, strong enough for washing away clot and foreign bodies from wounds, yet soft enough for washing the eyes. Ideal for the industrial sector, production industry, laboratories, transport of dangerous goods, first responders, law enforcement, firefighters, military and first rescue services, such as disaster relief. Plum Wound and Eyewash Spray is CE marked according to the Guidelines for Classification of Medical Device Directive (MDD93/42/EEC).

The product is available in the following sizes:

- 50 ml bottle (art. no. 45531). Rinse time: approx. 5 minutes. Carton: 24 x 50 ml
- 250 ml bottle (art. no. 4554). Rinse time: approx. 15 minutes. Carton: 12 x 250 ml



QuickRinse® eyewash ampoules

1 x 5 eyewash ampoules with 20 ml 0.9 % sterile sodium chloride. The ampoules provides easy, quick and hygienic relief for rinsing dust and dirt from the eyes. Suitable for refilling QuickSafe®, QuickFix&Rinse® and QuickRinse® stations.

The product is available in the following sizes:

- 5 x 20 ml ampoules (art. no. 5160). Carton: 4 x 5 ampoules.



Plum ON-the-GO Eyewash set

Practical eyewash set with 3 x 500 ml bottle containing 0.9% sterile sodium chloride solution. The Plum first aid eye wash is suitable for rinsing dirt (such as dust, dirt, metal and wood splinters, etc.) and even certain chemicals such as oils and solvents. It comes in an easy-to-use bottle, with an ergonomic cap, a dust cover and clear instructions for use. The box is easy to carry, you can take it with you in the car while you get to the doctor. The bottles can be used individually, but are also suitable for refilling eyewash stations. Shelf life: 3.5 years. Rinse time: approx. 15 minutes.

- Article number: 4620
- Packaging: 3 x 500 ml
- Rinse time: approx. 15 minutes



MUSEUM HILL CENTRE, MUSEUM ROAD
P.O. BOX 75534 00200
NAIROBI, KENYA
Tel: 0720 714 337
Email: info@medsurgehealth.co.ke
Web: www.medsurgehealth.co.ke