

NON-CONTACT TONOMETER (I-CARE)











Non-Contact Tonometer (I-Care)

Description

The first rebound tonometer introduced to the world by iCare. The original iCare tonometer is the right choice for a reliable basic device. The ease of use makes it an important instrument for general practitioners, optometrists, occupational health care, pharmacies and other medical personnel. It is handy in clinics for both adult and pediatric readings as well as in the field, i.e. in mass glaucoma screenings.







Key features

- Suitable for all patients
- Quick to measure
- Easy to use
- Durable
- Consistent and repeatable
- No drops, air or calibration needed

Specifications

	1 1		m	\sim	-	CI	_	-	0	
•	. ,	п		_		~ I	()		•	

Weight

Battery

Measurement range

Accuracy

 Repeatability (coefficient of variation) 13-32 mm x 45-80 mm x 230 mm

155 g (without batteries), 250 g (4 x AA batteries)

4 x AA non-rechargeable batteries

7 - 50 mmHg

(95 % tolerance interval relative to manometry):

±1.2 mmHg (20 mmHg) and ±2.2 mmHg (>20 mmHg)

<8%





MUSEUM HILL CENTRE, MUSEUM ROAD P.O. BOX 75534 00200 NAIROBI, KENYA Tel: 0720 714 337

Tel: 0720 714 337 Email: info@medsurgehealth.co.ke Web: www.medsurgehealth.co.ke

