

OFF-ROAD HEAVY DUTY POWER WHEELCHAIR

medsurge 
technics for life



OFF-ROAD HEAVY DUTY POWER WHEELCHAIR

DESCRIPTION

Using gaming-type cushion and backrest, the appearance of the line is textured, you can adjust the sitting and lying, and sitting for a long time is not tired



FEATURES

- Imported brand controller 60A,
- Taiwan Gear worm Motor.
- Electromagnetic brake motor , sensitive control, parking immediately be safer.
- Max weight load: 180kgs
- Optional used with captain seat.

Specifications

• Recline Available	Manual Reclining backrest
• Backrest Reclining	Max 135 degree
• Motors	DC500W*2 pcs
• Brake	Electromagnetic Brake
• Controller	60A
• Battery	24V 55AH Colloidal lead acid battery
• Charger	DC24V6A AC115-230V
• Charging time	6-8H
• Max speed	13km/h
• Driving distance	40km
• Product size	115*64*118cm
• Seat depth	48cm
• Seat width	50cm
• Seat height	62cm
• Front wheel	10" 3"
• Rear Wheel	14" 3"
• Safe load	180kg
• Climbing	8°-12°
• Package	98*70*84cm
• N.W / G.W	117kg /127kgs



MEDSURGE HEALTHCARE LIMITED

MUSEUM HILL CENTRE, MUSEUM ROAD

P.O. BOX 75534 00200

NAIROBI, KENYA

Tel: 0720 714 337

Email: info@medsurgehealth.co.ke

Web: www.medsurgehealth.co.ke