

AMBULANCE FOLDING STRETCHER

medsurge 
technics for life



AMBULANCE FOLDING STRETCHER

YDC-3K03 Ambulance Folding Stretcher

DESCRIPTION

Used in hospital or emergency center to transport patient.



FEATURES

- The main material is high-strength aluminum alloy.
- The height of whole stretcher can be adjusted in more than 4 positions, is convenient for CPR operation.
The height of front part and back part can be adjusted in 6 positions.
- The stretcher's leg folding mechanism and height adjustment is controlled through the front and rear pull rod operation.
- The guardrail using polymer plastics, the adjustable structure is easy for patient to up and down the stretcher, this stretcher adopts shoulder belts.
- It also equipped with a three-stage IV stand, the highest length is 80cm.
- The mattress is made of wear-resistant mesh cloth, concave shape, with a thickness of 10cm and disposable welding, flame retardant and waterproof.
- The universal wheels on the front end adopts a elastic latch to lock.
They can be locked to go straight or universal as needed.
- It's better to use F07 one-piece fixed device, the fixing device uses 3-point fixation, which is stable and secure.



Fixing device

Technical data

- Loading height: 56cm
- Wheel size: Φ 200mm
- Unfolded Size: 197 x 57 x 90cm
- High position 197 x 57 x 39cm
- Low position
- Packing Quantity: 1
- Packing Size: 204 x 62 x 52cm



MEDSURGE HEALTHCARE LIMITED
MUSEUM HILL CENTRE, MUSEUM ROAD
P.O. BOX 75534 00200
NAIROBI, KENYA
Tel: 0720 714 337
Email: info@medsurgehealth.co.ke
Web: www.medsurgehealth.co.ke