

Portable Pedal Manual Suction

medsurge 
technics for life



Portable Pedal Manual Suction

DESCRIPTION

The pedal suction unit is applicable for medical care unit as clinical suction under negative pressure.

FEATURES

- Product characteristics: small size, light weight, easy to operate, no power is required.
- Tread the stepping rod to drive the piston to move repeatedly to generate negative pressure. Use the non-return valve to generate positive pressure.
- The pedal suction unit is safe and reliable.
- Compact structure makes it easy to operate and convenient to carry.
- Medical suction with high negative pressure and high flow rate.

TECHNICAL SPECIFICATIONS

PARAMETER

SPECIFICATION

Ultimate Negative Pressure Valve	$\geq 0.08\text{MPa}$
Maximum Free-air Flow	$\geq 20\text{L/min}$
Reservoir Volume	$\geq 1000\text{ml}$, 1pc
Pedal Force	$\leq 245\text{N}$



MUSEUM HILL CENTRE, MUSEUM ROAD
P.O. BOX 75534 00200
NAIROBI, KENYA
Tel: 0720 714 337
Email: info@medsurgehealth.co.ke
Web: www.medsurgehealth.co.ke